

## HEALTH SURVEY SAMPLING INFORMATION

**1) General:** The Health Survey sampling kit consists of (4) four containers:

- (2) 1 quart sized plastic bottles for lead and inorganics analyses
- (1) 1 pint sized plastic bottle for nitrate and nitrite analysis
- (1) Sterile whirlpak container for bacteria analysis

**2) Sampling Container Preparation:**

Write name and sampling date on bottle labels and white area on sterilized whirlpak baggie. Use ink pen, being careful not to puncture the baggie.

**3) Site Selection:**

Sample sites/water spigots should be chosen carefully to obtain a representative sample of well water. However, precautions must be taken to choose a spigot that will not contaminate the sample. Spigots to be **avoided** include:

- 1) Spigots that are not used on a regular basis, with possible visible algae growth.
- 2) Spigots with leaking packing nuts or cracks.
- 3) Spigots with non-removable aerators, screens, or filters.
- 4) Spigots connected to plumbing in or around bathrooms.
- 5) Spigots that come **after** a whole-house softener or filter system.
- 6) Any hose. You can use a hose to flush the water, but remove it and collect directly from the spigot.
- 7) Pull-down kitchen faucets (sprayer and faucet in one).

The kitchen sink spigot is usually the best choice. However, for single lever faucets make sure **cold water only is run**. All aerators, screens, and filters should be removed or placed in bypass prior to sampling. New wells, or wells that have not been used recently, should be run before sampling to flush stagnant/contaminated water from the system.

**4) Sampling Procedure:**

One (1) quart plastic bottle for lead analysis needs to be filled first by **first draw** without letting the water flush too long. The EPA recommends a minimum stagnation time of 6 hours prior to sampling for lead. Ackuritlabs recommends sampling in the early morning or early evenings after work to ensure that stagnant water conditions exist. After sampling for lead (filling 1 quart container), let water continue to run for about five minutes. The water velocity should be reduced to a pencil sized stream prior to sampling and adjusted to avoid water from running around the threaded portion of the spigot head. Once the water velocity has been established, fill the other 1 quart and 1 pint bottle. The whirlpak sample baggie should be opened by tearing off the perforated top and using the paper tabs to open the baggie with care taken to avoid contact with the baggie lip. The whirlpak should then be filled without touching it to the spigot and then closed with a twirling motion flipped end over end using the weight of the water and secured shut using the wire tabs by twisting them together. The filled containers should be kept upright and cool until delivery to the lab.

**Bacteriological samples must be delivered to the lab within 24 hours of sample collection. Bacteriological samples are accepted Monday through Thursday and Friday only by prior arrangement and for an additional fee.**